Software Installation Guide:

BioTuring Colab

Introduction

- Software to tackle biomedical challenges

Colaboratory, or "Colab" for short, is a product from Bioturing. Colab has a variety of features and pre-built notebooks their user can download and use. We are providing many tools that help users to post their data and analyze the reports. Our product can be used to write and execute arbitrary python, R code, Golang, Julia, RStudio, VS Code and many more through the browser, and is especially well suited to data analysis and education. More technically, Colab is a hosted Jupyter notebook service that requires no setup to use. Users can build their own notebook.

System Requirements

Before installing the Colab, some pre-installation steps are required:

Basic Requirements for BioColab installation

	Basic recommendation	Optional				
		This is basic requirement to start BioColab				
		and based on requirement Resources as well				
CPU	16 core	as machine can be added				
RAM	64 Gb	As above				
HDD	/ partition can be 100 GB	as above				
	Data Volume : 1TB	As above				
	Any OS. Ubuntu 20.04 and above.	BioColab is more supportive with Linux OS. For better performance linux OS is recommended.				
Instance	Support any type of instance type. Depend on need	AWS g5xlarge in case using GPU				
Platform	Docker / Kubernetes					
Note						
	In case we want to use Notebook with BioColab. We need to install NVIDIA and below would be a requirement. The system has one or multiple NVIDIA GPU(s) (at least 16 GB memory per GPU) - with Turing architecture or above.					
	BioColan supports any Linux OS. We are recommending Ubuntu 20.04 or above.					
	SSL can be configured later also.					
	Please contact support@bioturing.com to get the token for your company.					
Security						
	The BioColab platform uses HTTPS protocol to securely communicate over the network.					
	All of the users need to authenticate using a BioTuring account or the company's SSO to access the platform.					
	We highly recommend setting up a private VPC network for IP restriction.					
	The data stays behind the company firewall.					
	The BioColab platform does not track any usage logs.					
Data visibility						

Data can be uploaded to Personal Workspace or Data Sharing group.		
In the Personal Workspace, only the owner can see and manipulate the data she/he uploaded.		
In the Data Sharing group, only people in the group can see the data.		
In the Data Sharing group, only people with sufficient permissions can manipulate the data.		

- The system has one or multiple NVIDIA GPU(s) (at least 16 GB memory per GPU) with Colab Bioturing architecture or above.
- The system is running Ubuntu 20.04 or above.
- SSL certificate and a domain name for users to securely access the platform on the web browser. I can be installed later too.
- Please contact support@bioturing.com to get the token for your company.

Note: The ideal system that we recommend for most companies is AWS g5.8xlarge for GPU based. Instance can be chosen based on requirement. If the notebook is not based on GPU, we can select lower.

Our Product is containerized applications, Here we are illustrated GPU based instances, Kindly select based on your requirement. It can be run on Docker using Docker engine and Kubernetes.

Download and Install

Note: We suggest starting from scratch to avoid package/driver conflicts.


For Tag naming conversion, kindly select based on your architecture.

Login to AWS console with admin user account to launch an Ec2 instance.

Note: It's up to the client. How they are going to manage infrastructure, Load, Network,

Access and traffic ...etc.

- Create a VPC
 - https://docs.aws.amazon.com/vpc/latest/userguide/what-is-amazon-vpc.h
 tml
 - https://docs.aws.amazon.com/vpc/latest/userguide/create-vpc.html
- 1. Select the appropriate region for VPC.
- Search VPC on the search box.
- 3. Click on Your VPCs New



4. Click on Create VPC





5. Follow the steps given in the image below.

VPC > Your VPCs > Create VPC

Create VPC Info


A VPC is an isolated portion of the AWS Cloud populated by AWS objects, such as Amazon EC2 instances.






Create Subnet









Verify Router Table










• IGW - Attach to VPC




• Update Router













• Create instance based on your requirement






Create a new key for SSH login to server






Select the VPC, which we created earlier.




• Create Security Group based on requirement


vpc-00c3c169 172.31.0.0/16





· Attach Disk for data volume





• Assign Elastic IP to this instance













Cancel

Associate




Associate Elastic IP address





• Verify Public IP of instance



Login to instance



Update the system:

- o sudo apt update && sudo apt upgrade -y
- sudo apt install build-essential wget curl gnupg lsb-release ca-certificates xfsprogs -y
- Install NVIDIA CUDA Toolkit 11.7.

- wget
 https://developer.download.nvidia.com/compute/cuda/11.7.1/local_instal
 lers/cuda_11.7.1_515.65.01_linux.run
- sudo sh cuda_11.7.1_515.65.01_linux.run

```
ubuntu@ip-192-168-0-199:~$ wget https://developer.download.nvidia.com/compute/cuda/11.7.1/local_installers/cuda_11.7.1_515.65.01_linux.run nux.ru--2023-03-30 08:20:01-- https://developer.download.nvidia.com/compute/cuda/11.7.1/local_installers/cuda_11.7.1_515.65.01_linux.run Resolving developer.download.nvidia.com (developer.download.nvidia.com)... 152.195.19.142 Connecting to developer.download.nvidia.com (developer.download.nvidia.com)|152.195.19.142|:443... nconnected. HTTP request sent, awaiting response... 200 OK Length: 3524358811 (3.3G) [application/octet-stream] Saving to: 'cuda_11.7.1_515.65.01_linux.run'
                                                    100%[============
cuda_11.7.1_515.65.01_linux.run
                                                                                                                                        2023-03-30 08:20:29 (119 MB/s) - 'cuda_11.7.1_515.65.01_linux.run' saved [3524358811/3524358811]
  buntu@ip-192-168-0-199:~$ sudo sh cuda_11.7.1_515.65.01_linux.run
ubuntu@ip-192-168-0-199:~$ sudo sh cuda_11.7.1_515.65.01_linux.run
Failed to verify gcc version. See log at /var/log/cuda-installer.log for details.ubuntu@ip-192-168-0-199:~$ cat /var/log/cuda-installer.log
[INFO]: Driver not installed.
[INFO]: Checking compiler version...
 [INFO]: gcc location:
[ERROR]: Missing gcc. gcc is required to continue.
ubuntu@ip-192-168-0-199:∼$ sudo apt install build-essential wget curl gnupg lsb-release ca-certificates xfsprogs -y
Reading package lists... Done
Building dependency tree
Reading state information... Done
lsb-release is already the newest version (11.1.0ubuntu2).
```

- Reboot the server and try again.
 - Ispci | grep -i nvidia
 - uname -m && cat /etc/*release

```
Last login: Thu Mar 30 08:17:02 2023 from 222.253.44.242

ubuntu@ip-192-168-0-199:~$ gcc --version
gcc (Ubuntu 9.4.0-1ubuntu1~20.04.1) 9.4.0

Copyright (C) 2019 Free Software Foundation, Inc.
This is free software; see the source for copying conditions. There is NO warranty; not even for MERCHANTABILITY or FITNESS FOR A PARTICULAR PURPOSE.

ubuntu@ip-192-168-0-199:~$ sudo sh cuda_11.7.1_515.65.01_linux.run
```

Accept license

End User License Agreement
----NVIDIA Software License Agreement and CUDA Supplement to
Software License Agreement. Last updated: October 8, 2021

The CUDA Toolkit End User License Agreement applies to the NVIDIA CUDA Toolkit, the NVIDIA CUDA Samples, the NVIDIA Display Driver, NVIDIA Nsight tools (Visual Studio Edition), and the associated documentation on CUDA APIs, programming model and development tools. If you do not agree with the terms and conditions of the license agreement, then do not download or use the software.

Last updated: October 8, 2021.

Preface

Do you accept the above EULA? (accept/decline/quit):

CUDA Installer
- [X] Driver
 [X] 515.65.01
+ [X] CUDA Toolkit 11.7
 [X] CUDA Demo Suite 11.7
 [X] CUDA Documentation 11.7
- [] Kernel Objects
 [] nvidia-fs
Options
Install

Up/Down: Move | Left/Right: Expand | 'Enter': Select | 'A': Advanced options

```
ubuntu@ip-192-168-0-199:~$ sudo sh cuda_11.7.1_515.65.01_linux.run
===========

Summary =
==========

Driver: Installed
Toolkit: Installed in /usr/local/cuda-11.7/

Please make sure that
- PATH includes /usr/local/cuda-11.7/bin
- LD_LIBRARY_PATH includes /usr/local/cuda-11.7/lib64, or, add /usr/local/cuda-11.7/lib64 to /etc/ld.so.conf and run ldconfig as root
To uninstall the CUDA Toolkit, run cuda-uninstaller in /usr/local/cuda-11.7/bin
To uninstall the NVIDIA Driver, run nvidia-uninstall
Logfile is /var/log/cuda-installer.log
ubuntu@ip-192-168-0-199:~$ ■
```

```
ubuntu@ip-192-168-0-199:~$ lspci | grep -i nvidia
00:1e.0 3D controller: NVIDIA Corporation Device 2237 (rev a1)
ubuntu@ip-192-168-0-199:~$ uname -m && cat /etc/*release
x86 64
DISTRIB_ID=Ubuntu
DISTRIB_RELEASE=20.04
DISTRIB_CODENAME=focal
DISTRIB_DESCRIPTION="Ubuntu 20.04.6 LTS"
NAME="Ubuntu"
VERSION="20.04.6 LTS (Focal Fossa)"
ID=ubuntu
ID_LIKE=debian
PRETTY_NAME="Ubuntu 20.04.6 LTS"
VERSION ID="20.04"
HOME_URL="https://www.ubuntu.com/"
SUPPORT URL="https://help.ubuntu.com/"
BUG REPORT URL="https://bugs.launchpad.net/ubuntu/"
PRIVACY POLICY URL="https://www.ubuntu.com/legal/terms-and-policies/privacy-policy"
VERSION CODENAME=focal
UBUNTU CODENAME=focal
ubuntu@ip-192-168-0-199:~$
```

nvidia-smi

ubuntu@ip-192-168-0-199:~ \$ nvidia-smi Thu Mar 30 08:36:04 2023								
NVIDIA-SMI 515.65.01	· ·							
GPU Name Persistence-M Bus-Id Disp.A Fan Temp Perf Pwr:Usage/Cap Memory-Usage	Volatile Uncorr. ECC							
0 NVIDIA A10G Off 00000000:00:1E.0 Off 00% 22C P0 57W / 300W 0MiB / 23028MiB	N/A							
+								

```
ubuntu@ip-192-168-0-199:/$ nvidia-smi
Thu Mar 30 10:11:57 2023
 NVIDIA-SMI 515.65.01
                          Driver Version: 515.65.01
                   Persistence-M
                                 Bus-Id
                                                Disp.A
                                                         Volatile Uncorr. ECC
                                                         GPU-Util Compute M.
 Fan Temp Perf Pwr:Usage/Cap
                                          Memory-Usage
                                                                       MIG M.
      NVIDIA A10G
                                  00000000:00:1E.0 Off
                     56W / 300W
       25C
                                    2992MiB / 23028MiB
               P0
                                                                      Default
                                                                          N/A
  GPU
        GΙ
             CI
                        PID
                              Type
                                     Process name
                                                                   GPU Memory
         ID
                       7829
    0
        N/A N/A
                                     ...t2d temp/tmpofbrlq31/exec
                                                                       249MiB
    0
        N/A N/A
                       7830
                                     ...t2d_temp/tmpofbrlq31/exec
                                                                       249MiB
        N/A N/A
                                 С
                                     ...t2d temp/tmpofbrlq31/exec
                                                                       249MiB
    0
                       7831
    0
        N/A N/A
                       7832
                                 C
                                     ...t2d_temp/tmpofbrlq31/exec
                                                                       249MiB
    0
        N/A N/A
                       8272
                                 С
                                     ...t2d_temp/tmp754zzxnp/exec
                                                                       249MiB
                                 C
    0
        N/A N/A
                                     ...t2d_temp/tmp754zzxnp/exec
                                                                       249MiB
                       8273
                                     ...t2d_temp/tmp754zzxnp/exec
                       8274
                                 C
    0
        N/A N/A
                                                                       249MiB
    0
        N/A N/A
                       8275
                                С
                                     ...t2d_temp/tmp754zzxnp/exec
                                                                       249MiB
                                     ...t2d_temp/tmp754zzxnp/exec
    0
        N/A N/A
                       8276
                                                                       249MiB
    0
        N/A N/A
                                     ...t2d_temp/tmp754zzxnp/exec
                                                                       249MiB
                       8278
    0
        N/A N/A
                       8280
                                 С
                                     ...t2d temp/tmp754zzxnp/exec
                                                                       249MiB
        N/A N/A
                       8281
                                     ...t2d_temp/tmp754zzxnp/exec
                                                                       249MiB
ubuntu@ip-192-168-0-199:/$
```

- Install docker.
 - curl https://get.docker.com | sh
 - sudo systemctl --now enable docker

Method 1: Using Docker

• Configure DNS entry. We are using Cloudflare.



Install the NVIDIA container toolkit.

o distribution=\$(./etc/os-release;echo \$ID\$VERSION_ID) \
&& curl -fsSL https://nvidia.github.io/libnvidia-container/gpgkey | sudo
gpg --dearmor -o /usr/share/keyrings/nvidia-container-toolkit-keyring.gpg
\
&& curl -s -L
https://nvidia.github.io/libnvidia-container/\$distribution/libnvidia-containe
r.list | \
sed 's#deb https://#deb
[signed-by=/usr/share/keyrings/nvidia-container-toolkit-keyring.gpg]
https://#g' | \

- o sudo tee /etc/apt/sources.list.d/nvidia-container-toolkit.list
- sudo apt update
- sudo apt install nvidia-docker2
- sudo systemctl restart docker
- nvidia-docker version

```
ubuntu@ip-192-168-0-199:-$ nvidia-docker version

NVIDIA Docker: 2.12.0

Client: Docker Engine - Community

Version: 23.0.2

API version: 1.42

60 version: g01.19.7

Git commit: 569dd73

Built: Mon Mar 27 16:16:18 2023

OS/Arch: linux/amd64

Context: default

permission denied while trying to connect to the Docker daemon socket at unix:///var/run/docker.sock: Get "http://%2Fvar%2Frun%2Fdocker.sock/v1.24/version": dial unix /var/run/docker.sock: Get "http://%2Fvar%2Frun%2Fdocker.sock/v1.24/version": dial unix /var/run/docker.sock: Onnect: permission denied ubuntu@ip-192-168-0-199:-$ 

□
```

Issue:

permission denied while trying to connect to the Docker daemon socket at unix:///var/run/docker.sock: Get "http://%2Fvar%2Frun%2Fdocker.sock/v1.24/version": dial unix /var/run/docker.sock: connect: permission denied ubuntu@ip-192-168-0-199:~\$ Is -lhart /var/run/docker.sock srw-rw---- 1 root docker 0 Mar 30 08:59 /var/run/docker.sock

```
permission denied while trying to Connect to the bocker daemon unix /var/run/docker.sock: connect: permission denied bubuntu@ip-192-168-0-199:-$ 1s -lhart /var/run/docker.sock sww.nv---1 root docker 0 Mar 30 08:59 /var/run/docker.sock ubuntu@ip-192-168-0-199:-$ sudo chmod 666 /var/run/docker.sock ubuntu@ip-192-168-0-199:-$ nvidia-docker version
ubuntuelp-192-188-0-199:\( \infty\) nviola-
NVIDIA Docker: 2.12.0

Client: Docker Engine - Community
Version: 23.0.2

API version: 1.42

Go version: gol.19.7

Git commit: 569dd73
                                              Mon Mar 27 16:16:18 2023 linux/amd64
 Context:
                                              default
 Server: Docker Engine - Community
                                             23.0.2
1.42 (minimum version 1.12)
go1.19.7
219f21b
   API version:
Go version:
Git commit:
                                              Mon Mar 27 16:16:18 2023 linux/amd64
   Built:
  Experimental:
containerd:
Version:
                                              false
                                             1.6.19
1e1ea6e986c6c86565bc33d52e34b81b3e2bc71f
   GitCommit:
                                              1.1.4
v1.1.4-0-g5fd4c4d
    Version:
 GitCommit:
docker-init:
   Version:
GitCommit:
                                              0.19.0
```

- Make sure that /dev/shm size is at least half of physical memory.
 To change the configuration for /dev/shm, add one line to /etc/fstab. For example, if the system has 128 GB of physical memory:
 - tmpfs /dev/shm tmpfs defaults,size=64g 0 0
- Run the command below to make the change immediately:
 - sudo mount -o remount /dev/shm

```
ubuntu@ip-192-168-0-199:~$ sudo vi /etc/fstab
ubuntu@ip-192-168-0-199:~$ sudo mount -o remount /dev/shm
```

```
ubuntu@ip-192-168-0-199:~$ df -h
Filesystem
                 Size
                        Used Avail Use% Mounted on
/dev/root
                 146G
                         13G
                              133G
                                     9% /
devtmpfs
                  63G
                           0
                               63G
                                     0% /dev
tmpfs
                  64G
                               64G
                                     0% /dev/shm
                           0
tmpfs
                  13G
                        1.1M
                               13G
                                     1% /run
tmpfs
                                     0% /run/lock
                 5.0M
                           0
                              5.0M
tmpfs
                  63G
                           0
                               63G
                                     0% /sys/fs/cgroup
/dev/loop0
                  25M
                         25M
                                 0 100% /snap/amazon-ssm-agent/6312
/dev/loop1
                                 0 100% /snap/canonical-livepatch/164
                 9.0M
                        9.0M
/dev/loop2
                                 0 100% /snap/core20/1852
                  64M
                         64M
/dev/nvme0n1p15
                                     6% /boot/efi
                 105M
                        6.1M
                               99M
/dev/loop3
                 117M
                        117M
                                 0 100% /snap/core/14946
/dev/loop4
                                 0 100% /snap/lxd/24061
                  92M
                         92M
/dev/loop5
                                 0 100% /snap/snapd/18596
                  50M
                         50M
                                 0 100% /snap/core18/2714
/dev/loop6
                  56M
                         56M
tmpfs
                  13G
                           0
                                     0% /run/user/1000
                               13G
/dev/loop7
                  44M
                         44M
                                 0 100% /snap/certbot/2836
ubuntu@ip-192-168-0-199:~$ 🗍
```


- o sudo mkdir -p /config/ssl
- sudo mv tls.crt /config/ssl
- sudo mv tls.key /config/ssl
- Create default directories to store user data. We highly recommend using persistent storage for these directories. In the commands below, we use an empty EBS volume.
- Pull the BBrowserX image.
 - sudo mkfs -t ext4 /dev/nvme2n1
 - sudo mkdir /data
 - o sudo mount /dev/nvme2n1 /data
 - sudo mkdir /data/app_data
 - sudo mkdir /data/user_data

sudo docker pull bioturing/bioturing-colab:1.0.1

• Run the docker image.

```
docker run -it -d \
-e WEB_DOMAIN='<yourcompany.com>' \
-e BIOTURING_TOKEN='<your token from BioTuring>' \
-e SSO_DOMAINS='<your company email address, example: @bioturing.com>' \
-e ADMIN_USERNAME='<your admin user name>' \
-e ADMIN_PASSWORD='<your admin password>' \
-v /data/user_data:/data/user_data \
-v /data/app_data:/data/app_data \
-name bioturing-colab \
--gpus all \
--shm-size=64gb \
-p 443:443 \
-p 80:80 \
bioturing/bioturing-colab:1.0.1
```

Wait for a few minutes for the platform to download all of the required services. After that, the BioTuring Colab is up and running.



Method 2: Using Kubernetes

====== Reference

https://www.howtoforge.com/how-to-install-containerd-container-runtime-on-ubuntu-22-04/#comments

https://blog.antosubash.com/posts/setup-micro-k8s-with-ubuntu

• Patch container engines (Docker, Containerd)

Install NVidia container toolkit on each node following the guide: https://docs.nvidia.com/datacenter/cloud-native/container-toolkit/install-guide.html

Check container engines (Docker, Containerd)

https://github.com/bioturing/installation

• Check container engines (Docker, Containerd)

For microk8s:

o microk8s kubectl describe no | grep Runtime

For vanilla:

- kubectl describe no | grep Runtime
- If container engine is Containerd, add these lines to:
 /etc/containerd/config.toml

```
privileged_without_host_devices = false
base_runtime_spec = ""
[plugins."io.containerd.grpc.v1.cri".containerd.runtimes.runc.options]
  SystemdCgroup = true
[plugins."io.containerd.grpc.v1.cri".containerd.runtimes.nvidia]
  privileged_without_host_devices = false
  runtime_engine = ""
  runtime root = ""
  runtime_type = "io.containerd.runc.v1"
  [plugins."io.containerd.grpc.v1.cri".containerd.runtimes.nvidia.options]
  BinaryName = "/usr/bin/nvidia-container-runtime"
  SystemdCgroup = true
[plugins."io.containerd.grpc.v1.cri".cni]
bin_dir = "/opt/cni/bin"
conf_dir = "/etc/cni/net.d"
```

After that, restart containerd

sudo systemctl restart containerd

- o sudo nvidia-container-cli --load-kmods info
- If container engine is Docker, add these lines to : /etc/docker/daemon.json

```
"default-runtime": "nvidia",

"runtimes": {
    "nvidia": {
        "path": "nvidia-container-runtime",
        "runtimeArgs": []
     }
}
```

- After that, restart docker
 - o sudo systemctl restart docker
 - o sudo nvidia-container-cli --load-kmods info

```
<mark>ountu@ip-192-168-0-186:∼$ m</mark>icrok8s kubectl describe no | grep Runtime
  Container Runtime Version: containerd://1.6.8
ubuntu@ip-192-168-0-186:~$ cat /etc/docker/daemon.json
    "runtimes": {
    "nvidia": {
            "path": "nvidia-container-runtime",
            "runtimeArgs": []
ubuntu@ip-192-168-0-186:~$ docker -v
Docker version 23.0.2, build 569dd73
ubuntu@ip-192-168-0-186:~$ sudo systemctl restart docker
do nvidia-container-cli --load-kmods infoubuntu@ip-192-168-0-186:~$ sudo nvidia-container-cli --load-kmods info
NVRM version: 515.65.01
CUDA version: 11.7
Device Index:
Device Minor:
Model:
                NVIDIA A10G
Brand:
                Nvidia
GPU UUID:
                GPU-fd9ffd32-13f4-61dd-ce58-197c9f1b9fb5
Bus Location:
                00000000:00:1e.0
Architecture:
```

Install BioTuring Colab on K8S

We support all k8s engines: GKE (Google Kubernetes Engine), EKS (Amazon Elastic Kubernetes Service), AKS (Azure Kubernetes Service), MicroK8s, and vanilla K8S.

Ensure that helm (version 3) is installed.

First, check the Helm version

- o microk8s enable helm3
- o microk8s helm3 version

Add BioTuring Helm charts

https://bioturing.github.io/charts/: for kubernetes

Example:

For Vanilla K8s:

- helm repo add bioturing https://bioturing.github.io/charts for Microk8s:
- microk8s helm3 repo add bioturing https://bioturing.github.io/charts

```
ubuntu@ip-192-168-0-186:~$ microk8s enable helm3

Infer repository core for addon helm3

Addon core/helm3 is already enabled
ubuntu@ip-192-168-0-186:~$ microk8s helm3 version
version.BuildInfo{Version:"v3.9.1+unreleased", GitCommit:"7112315b8d78eab23e9542c4fd824375429ca965", GitTreeState:"clean", GoVersion:"go1.19.5"}
ubuntu@ip-192-168-0-186:~$ microk8s helm3 repo add bioturing https://bioturing.github.io/charts/apps/
"bioturing" already exists with the same configuration, skipping
ubuntu@ip-192-168-0-186:~$

■
```

Helm chart Values Kubernetes: >=

Key	Туре	Default	Description
image.tag	string	"1.0.1"	image tag
secret.data.domain	string	"colab.com"	your domain
secret.data.ssodomains	string	1111	allow domains
secret.data.bbtoken	string	""	bioturing access token
secret.data.allowips	string	""	allow ips
secret.admin.username	string	admin	username
secret.admin.password	string	turing2022	password
secret.server.useletsencrypt	string	"false"	
secret.server.lcall	string	"C.UTF-8"	
secret.server.lclang	string	"C.UTF-8"	
secret.server.certificate	string	""	CRT base64 string
secret.server.key	string	""	KEY base64 string
service.type	string	ClusterIP	
service.ports.http.port	int	80	
service.ports.https.port	int	443	
persistence.dirs.app.size	string	5Gi	APP size
persistence.dirs.app.storageClass	string	1111	
persistence.dirs.user.size	string	5Gi	USER size
persistence.dirs.shm.size	string	1Gi	SHM size
persistence.dirs.user.storageClass	string	""	
persistence.dirs.user.existingClaim	bool	FALSE	
ingress.enabled	bool	TRUE	
ingress.className	string	1111	
ingress.annotations	object	{}	
ingress.tls.enabled	bool	TRUE	
resources	object	{}	
autoscaling	object	{}	

nodeSelector	object	{}	
tolerations	object	{}	
affinity	object	{}	
podAnnotations	object	{}	
podSecurityContext	object	{}	
securityContext	object	{}	
serviceAccount.name	string	""	
gpu.enabled	bool	TRUE	
gpu.runtimeClassName	string	"nvidia"	

• For Containerd runtime :

gpu.runtimeClassName="nvidia"

• For Docker runtime :

gpu.runtimeClassName=""

• Simple Installation (Recommended):

```
End User License Agreement

NVIDIA Software License Agreement and CUDA Supplement to
Software License Agreement. Last updated: October 8, 2021

The CUDA Toolkit End User License Agreement applies to the
NVIDIA CUDA Toolkit, the NVIDIA CUDA Samples, the NVIDIA
Display Driver, NVIDIA Nsight tools (Visual Studio Edition),
and the associated documentation on CUDA APIs, programming
model and development tools. If you do not agree with the
terms and conditions of the license agreement, then do not
download or use the software.

Last updated: October 8, 2021.

Preface

Do you accept the above EULA? (accept/decline/quit):
accept
```

```
CUDA Installer
- [X] Driver
[X] 515.65.01
+ [X] CUDA Toolkit 11.7
[X] CUDA Demo Suite 11.7
[X] CUDA Documentation 11.7
- [] Kernel Objects
[] nvidia-fs
Options
Install

Up/Down: Move | Left/Right: Expand | 'Enter': Select | 'A': Advanced options
```

```
ubuntu@ip-192-168-0-186:-$ bash ./install.k8s.sh
Your K8S engine [vanilla, microk8s]: microk8s
Do you need install CUDA Toolkit [y, n]: n
Ignore re-install CUDA
```

Going through this interactive installation to finish the installation. After this step, just access the BioTuring System via the specified domain in the installation process. If it's not in the DNS, please add the ip/domain to the local machine DNS host file.

Check pods information

- microk8s kubectl get all
- microk8s kubectl get pods
- microk8s kubectl get services --all-namespaces
- microk8s kubectl get services
- microk8s kubectl get pvc
- microk8s kubectl logs bioturing-colab
- microk8s.kubectl -n ingress get pods
- microk8s.kubectl -n ingress logs <your pod name here> | grep reload

Check secrets

- bioturing-colab-tls
- bioturing-colab
- bioturingregred

microk8s kubectl edit secrets mysecret

Example:

microk8s kubectl edit secrets bioturing-colab-tls

Manual Installation

Please replace paths to your certificate, key, admin password, and other helm chart values of your choice.

BBTOKEN="USE TOKEN OBTAINED FROM BIOTURING"

SSLCRT="base64 -w 0 ./bioturing.com.crt" # <- (REPLACE THIS WITH A PATH TO YOUR CRT CERTFICATE)

SSLKEY="base64 -w 0 ./bioturing.com.key" # <- (REPLACE THIS WITH A PATH TO YOUR KEY)

ADMIN USERNAME="admin"

ADMIN_PASSWORD="admin" # <- (CHANGE YOUR PASSWORD IF NECESSARY)

USELETSENCRYPT="false"

SVHOST="k8stest.bioturing.com" # <- (CHANGE THIS TO YOUR K8S INGRESS DOMAIN)

APP DATA SIZE="50Gi" # <- (CHANGE THIS TO YOUR APP-PVC SIZE)

USER_DATA_SIZE="100Gi" # <- (CHANGE THIS TO YOUR USER-PVC SIZE)

SHM SIZE="64Gi" # <- (CHANGE THIS TO YOUR SHM SIZE)

CHART VERSION="1.0.23" # <- (CHANGE IT IF NECESSARY)

LC_ALL="C.UTF-8" # <- (CHANGE IT IF NECESSARY)

LC_LANG="C.UTF-8" # <- (CHANGE IT IF NECESSARY)

For Microk8s:

microk8s helm3 repo update

microk8s helm3 registry login -u admin registry.bioturing.com

microk8s helm3 upgrade --install --set secret.data.bbtoken="\${BBTOKEN}" \

- --set secret.data.domain="\${SVHOST}" \
- --set secret.server.certificate="\${SSLCRT}" \
- --set secret.server.key="\${SSLKEY}" \
- --set secret.server.useletsencrypt="\${USELETSENCRYPT}" \
- --set secret.server.lcall="\${LC ALL}" \
- --set secret.server.lclang="\${LC LANG}" \
- --set secret.admin.username="\${ADMIN_USERNAME}" \

```
--set secret.admin.password="${ADMIN_PASSWORD}" \
--set persistence.dirs.app.size="${APP_DATA_SIZE}" \
--set persistence.dirs.user.size="${USER_DATA_SIZE}" \
--set persistence.dirs.shm.size="${SHM_SIZE}" \
bioturing bioturing/colab --version ${CHART_VERSION}
```

For Vanilla k8s:

```
helm repo update
helm registry login -u admin registry.bioturing.com
helm upgrade --install --set secret.data.bbtoken="${BBTOKEN}" \
--set secret.data.domain="${SVHOST}" \
--set secret.server.certificate="${SSLCRT}" \
--set secret.server.key="${SSLKEY}" \
--set secret.server.useletsencrypt="${USELETSENCRYPT}" \
--set secret.server.lcall="${LC_ALL}" \
--set secret.server.lclang="${LC_LANG}" \
--set secret.admin.username="${ADMIN_USERNAME}" \
--set secret.admin.password="${ADMIN_PASSWORD}" \
--set persistence.dirs.app.size="${APP_DATA_SIZE}" \
--set persistence.dirs.user.size="${SHM_SIZE}" \
--set persistence.dirs.shm.size="${SHM_SIZE}" \
bioturing bioturing/colab --version ${CHART_VERSION}
```

SSO setup: There are various service providers, who work as IDP. Here we use Jumpcloud and okta.

- please select yourself.

https://colab.bioturing.com/dashboard/sso

SSO configuration: kindly contact support@bioturing.com for all product related questions, issues, including training.

Configuration

All set with SSO

Troubleshooting

[Provide tips for troubleshooting common issues that may arise during or after installation.]

Issue:

permission denied while trying to connect to the Docker daemon socket at unix:///var/run/docker.sock: Get "http://%2Fvar%2Frun%2Fdocker.sock/v1.24/version": dial unix /var/run/docker.sock: connect: permission denied

Solution:

```
docker ps -a
groups
sudo groupadd docker
sudo groups
sudo usermod -aG docker $USER
sudo systemctl start docker
newgrp docker
docker ps -a

ubuntu@ip-192-168-0-199:~$ Is -lhart /var/run/docker.sock
srw-rw---- 1 root docker 0 Mar 30 08:59 /var/run/docker.sock
```

ubuntu@ip-192-168-0-199:~\$ sudo chmod 666 /var/run/docker.sock

Conclusion

